



NICON

Oxford Foot Model Training Day

3rd July 2008




Introducing...



The Oxford Foot Model

NICON



Why use a foot model?

- The foot contains 26 bones, 33 joints, 107 ligaments and 19 muscles.
- Intricate mechanics – particular issue where deformity present




NICON



Foot Anatomy and Function

NICON



Motivation for Development

- Original model proposed by Melissa Carson and Marian Harrington
- Measure foot deformity and outcome of treatment in Club Feet




NICON



Evolution of the OFM

- Aims:
 - Compatible with known foot function
 - Repeatable
 - Adaptable to foot deformity

NICON



Evolution of the OFM

- Include stance and swing phase
- Hallux definition
- Axis definitions
- Marker placement
- Static trial

VICON



Children and foot deformity

Challenges posed by:

- small surface area of the foot
- greater variability of foot motion
- deformities present in pathological feet



VICON



Published Studies


- Healthy Adults
Carson et al, *J Biomech* 2001
- Club Feet
Theologis et al, *JBJS* 2003
- Children
Stebbins et al, *Gait and Posture* 2006



VICON




Implemented as Vicon PlugIn




The Oxford Foot Model

- 3 segment model with optional hallux



Unique characteristics


- 4 markers per segment
- Not reliant on neutral position
- Not reliant on x-ray
- Primary axes on long axis of foot
- Compatible with PlugInGait
- Easily understood



VICON

Understanding the OFM

- Described primarily by planes
- Origin + 3 axes




VICON




Description of axes

- Tibia
 - Described by frontal plane




VICON



- Origin: AJC
- Primary axis: AJC – KJC
- Medio-lateral axis: Bimalleolar axis
- 3rd Axis: mutually perpendicular
- Same as for PlugIn Gait



- Hindfoot
 - Mid-sagittal plane of calcaneus




VICON



- Origin: HEEL marker
- Primary Axis: Parallel to floor and in plane of HEE, PCA and midpoint between STL and LCA
- Medio-Lateral Axis: Perpendicular to this plane
- 3rd Axis – mutually perpendicular


- Midfoot
 - Linking mechanism only



VICON

- Forefoot

- Transverse plane of metatarsals



VICON



- Origin: Midpoint between P5M and P1M
- Primary Axis: Origin to TOE projected into plane of P5M, D5M and D1M
- Vertical Axis – Perpendicular to this plane
- 3rd Axis – Mutually perpendicular

VICON



- Hallux
- Vector: D1M - HLX



VICON



Calculation of Angles

- Grood and Suntay* sequence
 - Flexion extension about transverse axis of proximal segment
 - Axial rotation about long axis of distal segment
 - Ab/Adduction about “floating” mutually perpendicular axis
 - NB – turn the corner at the ankle

* Grood ES, Suntay WJ. A joint coordinate system for the clinical description of three-dimensional motions: application to the knee. J Biomech Eng 1983;105:136-44.

VICON



Demonstration

- Foot clinical examination
- Marker placement
- Practice!



VICON



Demonstration

- Collecting and processing OFM

VICON



Practical Guidelines

- Camera set-up required



VICON



Practical Guidelines

- Who is suitable?
 - Age
 - Co-operation
 - Foot question
 - Degree of deformity
- Walking Aides
- Markers staying on

VICON



Practical Guidelines

- Processing
 - Tracking markers
 - Labelling
 - Checking for "sensibleness"
- Displaying Results

VICON



Practical Guidelines

- Interpreting Graphs
 - How many trials?
 - Consistency
 - What is abnormal?

VICON




Interpretation of Graphs

- Case Studies

VICON



Graphs



VICON

Case 1: Cavo-varus Foot (C769B)

- 9 year old boy with Cerebral Palsy
 - (R Hemiplegia) – considering surgery
- Independent community ambulator
- Uses R AFO and 1cm shoe raise
- Previously received B-toxin injections
- Walks up to 15 min, PE at school
- Walking recently deteriorated
- Pain in R ankle and heel

VICON



Clinical Examination


- (R) 5° hip flexion contracture
- (R) 50° popliteal angle
- (R) 10° equinus contracture with flexible cavo-varus foot posture (WB)
- 1.5cm leg length discrepancy (R shorter)

VICON




Lower Body Kinematics

- Type IV



VICON

Foot Model Kinematics



VICON

Conclusions

- Equinus from both hindfoot and forefoot
- Combination of both hindfoot varus AND forefoot supination (unusual)
- Adduction of foot from forefoot level
- Reduced ankle ROM
- Tibial torsion and hamstring tightness also contributing

VICON



Recommendations

- Option 1- Serial casting and orthotics, would need to correct flexible cavus and forefoot adduction
- Option 2- Surgery: gastroc lengthening, tib post lengthening, split tib ant transfer, derotation tibia, hamstring lengthening
- **Option 3- early surgery: gastroc and tib post lengthening**

VICON



Case 2: Planovalgus Foot (C303F)

- 17 year old boy with Cerebral Palsy/ Spastic Diplegia who walks with crutches or a Kaye walker and GRAFOs
- Previous multi-level soft tissue surgery in 2000
- Referred for consideration of further surgery

VICON




Clinical Examination

- 35° hip flexion contractures
- 5° (R) and 10° (L) knee flexion contractures
- 15° (R) and 25° (L) equinus contractures
- Uncorrectable planovalgus feet with mid-foot breaks


VICON



Lower Limb Kinematics



VICON



VICON



Conclusions

- (L) hindfoot plantarflexion and forefoot dorsiflexion
- Reduced ankle ROM
- Bilateral hindfoot eversion and forefoot supination
- Bilateral hindfoot internal rotation and forefoot abduction
- Reduced arch height (R)>(L)

VICON



Recommendations

- Subtalar arthrodesis and calcaneal lengthenings with gastroc lengthenings
 - + multi-level surgery

VICON



Case 3: Plano-valgus Foot (C521B)

- 10 yr old with (R) hemiplegia (CP)
 - but bilateral planovalgus feet
- Diagnosis unclear
- Wears bilateral AFOs
- No functional limitations

VICON




Clinical Examination

- (R) hip flexion contracture of 5°
- (R) popliteal angle 70°
- (R) achieve plantargrade position, with fixed planus and midfoot break, flexible valgus, hallux valgus
- Leg length discrepancy (2 cm shorter R)

VICON




Lower Body Kinematics



VICON

Foot Model Kinematics



VICON

Conclusions

- Dynamic foot deformity not as obvious as in static assessment
- Dynamic insufficiency at ankle level
- Compromised hallux function

VICON



Recommendations

- (R) side: AFO
- (L) side: AFO / foot orthosis
- may be surgical candidate if foot deformity deteriorates on either side.

VICON



Case 4: – Pre/Post Surgery

HISTORY

- 8 year old girl with R Hemiplegia
- Community ambulator
- Wears R fixed AFO and night splint
- Previous Botulinum Toxin injections and serial casting with some success


VICON



Clinical Examination – R side


- Hip flexion contracture - 15°
- Excessive anteversion - 25°
- Knee flexion contracture - 15°
- Fixed equinus - 20°
- Excessive external tibial torsion
- Mild Rectus and Gastrocnemius spasticity
- Poor strength and selective control at foot
- 3cm leg length discrepancy (R shorter)

Pre-op kinematic graphs



VICON

Oxford Foot Model



Surgery

- Plantarflexor lengthening
- Tibialis Posterior lengthening
- Split Tibialis Anterior transfer




VICON




Conventional graphs – outcome

“Bold is Old”



VICON

Foot model – outcome of surgery



Clinical Implications and Research

- Repeatability
 - Healthy
 - CP
- Findings from clinical practice
- Further work

VICON



Repeatability

- Intra/Inter-tester: Healthy Adults
Carson et al, *J Biomech* 2001
- 2 testers/ 2 subjects
- 1° to 6.5° variability

VICON



Repeatability


- Intra-tester: Healthy Children
Stebbins et al, *Gait and Posture* 2006
- 1 tester/ 15 subjects
- 1° to 7.5° variability

VICON



Repeatability

- Inter/intra-tester: healthy adult
- 5 testers
- 1 subject
- 0.5° - 7°




VICON



Repeatability

- Intra-subject: Children with CP



VICON




Clinical Findings

VICON



OFM Data ...

From 2004 to present: 163 patients



Cerebral Palsy - Diplegia

- 23 patients (43 feet) assessed



VICON



Cerebral Palsy - Diplegia

- Foot model information used to:
 - Determine cause of in-toeing/out-toeing
 - Determine level of foot drop
 - Guide orthotic intervention
 - Confirm visual impression of foot deformity
 - Aid clinical assessment

VICON



Cerebral Palsy - Diplegia

- Significant difference compared to healthy control data
 - ROM sagittal hindfoot
 - ROM sagittal forefoot
- Coronal and transverse planes:
 - Equal numbers in both directions
 - Diversity of foot deformity

VICON



Cerebral Palsy - Hemiplegia

- 70 patients (6 – 38 years)
- 6 bilateral (just started)



VICON



Cerebral Palsy - Hemiplegia


- Foot model information used to:
 - Assess dynamic foot motion
 - Assess outcomes of treatment
 - Monitor progression of foot deformity
 - Clarify controversial findings from lower limb kinematics
 - Determine level of foot drop
 - Corroborate clinical findings

VICON



Cerebral Palsy - Hemiplegia

- Significant difference with healthy controls
 - ↓ ROM hindfoot and forefoot sagittal
 - ↓ Maximum hindfoot and forefoot dorsiflexion
 - ↑ Hindfoot internal rotation



VICON




Club Foot

- 24 patients (7 female and 17 male, age range 6 to 24 years)
- 13 bilateral, 6 right and 5 left clubfeet = total of 37 feet
- Early age posterior-medial release

VICON



Results of Clubfoot Surgery



VICON



Club Foot

- The findings from the OFM were used to identify:
 1. Level of dynamic foot deformity (hindfoot and/or forefoot)
 2. Specify the type of surgery required
 3. Justify type of casting appropriate
 4. Clarify the source of foot rotation (tibia, hindfoot or forefoot)
 5. Corroborate clinical findings
- The OFM directly influenced the future management recommendations in 45% of patients seen

VICON



Other Populations

Neuro- 23 patients

- Polio
- Multiple Sclerosis
- HMSN
- Charcot Marie Tooth
- Nerve lesion
- Stroke
- Head Injury
- Spina Bifida

Ortho – 22 patients

- Rotational Mal-alignment
- In-toeing
- Flat feet
- Pain (hip, knees, feet)

- Toe walkers

VICON



Hemiplegic cerebral palsy: foot kinematics of 'good' foot

VICON



Method

- Subjects
 - 10 children (8-17 years)
 - Hemiplegic cerebral palsy
 - 6 pre-op, 4 post-op to affected side
 - Mildly affected (6 type I)
- Looked at
 - clinical exam
 - foot model (static & dynamic)

VICON




Clinical exam

- Non weight-bearing
 - Neutral (9 subjects)
- Weight-bearing
 - Planovalgus (5 subjects)
 - Valgus (2 subjects)


VICON



Static foot model



Dynamic foot model



Ongoing Development

- Hallux varus/valgus
- Normal static weight-bearing foot position compared to clinical exam and kinematics
- Quantifying skin/marker movement on the foot – use of markers
- Functional calculation of ankle joint complex axes
- Correlate with imaging

VICON



Questions?



VICON



Contributors to the Oxford Foot Model

- Dr Melissa Carson
- Dr Marian Harrington
- Dr Amy Zavatsky
- Mr Tim Theologis
- Mrs Nicky Thompson
- Dr Claudia Giacomozzi

VICON



- Dr Richard Baker
- Mr Timothy Pitt
- The Oxford Gait Lab team
- Vicon



VICON



VICON

Thank you

