

## MARKER PLACEMENT

**Patient Preparation** - see *6.4 Clinical Examination Protocol*

**Models** - standard marker placement protocol for:

- |                                       |         |   |
|---------------------------------------|---------|---|
| • Lower body Plug-in-Gait (PiG)       | 19 mkrs | 13 dynamic<br>4 wand (thigh x2, shank x2)<br>2 static |
| • Upper body Plug-in-Gait (UL PiG)    | 23 mrks | 23 dynamic  |
| • Oxford Trunk Model (OTM)            | 4 mkrs  | 3 dynamic<br>1 static (1 same as UL PiG)              |
| • Oxford Foot Model (OFM)             | 24 mkrs | 18 dynamic<br>2 wand (C-Peg x2)<br>4 static           |
| • Knee and hip functional calibration | 12 mkrs | 12 dynamic (6 same as OFM)                            |

**Combined marker sets:**

- |                             |         |  |
|-----------------------------|---------|--|
| • PiG, functional           | 31 mkrs | 25 dynamic<br>4 wand (thigh x2, shank x2)<br>2 static            |
| • PiG, functional, OTM      | 35 mkrs | 28 dynamic<br>4 wand (thigh x2, shank x2)<br>3 static            |
| • PiG, OFM, functional      | 49 mkrs | 37 dynamic<br>6 wands (thigh x2, shank x2, C-Peg x2)<br>6 static |
| • PiG, OFM, functional, OTM | 53 mkrs | 40 dynamic<br>6 wands (thigh x2, shank x2, C-Peg x2)<br>7 static |

Description of each model's marker placement is below, with images at the end of the document.

**Default models to use are:**

**Barefoot  
Shod**

**PiG, functional, OFM\*  
PiG, functional**

\* OFM may not be possible for small and/or clammy feet

The text below for PiG is adapted from Vicon documentation. Where left side markers only are listed; the positioning is identical for the right side. Forefoot and anterior tibial markers are placed in sitting. All other markers, where possible are placed with the patient standing.

### LOWER BODY PLUG-IN-GAIT (PiG)

Standard clinical marker set shown in bold; 19 markers total (2 removed for dynamic trials)

#### *Pelvis markers*

<b>LASI</b>	Left ASIS	Palpate from below anterior superior iliac spine. Place marker just above projection (lip) of ASIS.
<b>LISH</b>	Left iliac crest	Palpate along the iliac crest. Place marker at the midpoint between the SACR and LASI

In obese patients it can be difficult to place markers accurately on the ASIS and excessive soft tissue can obscure the markers. In this case, use the Hvidovre wand markers as a first option. If these will be knocked when walking, then move each ASIS marker laterally by an equal amount along the LASIS-RASIS axis. If it is not possible to move markers laterally without also moving them antero-posteriorly the ant-post ASIS to trochanter distance should be recorded and input into the model. The true inter-ASIS distance must then be recorded and entered on the PRF form.

<b>SACR</b>	Sacrum	Placed on the skin mid-way between the posterior superior iliac spines. Check positions of ASIS and SACR markers from the side, to make sure pelvic tilt represented by markers looks sensible.
-------------	--------	---

If losing visibility of the sacral marker is a problem, a skin marker can be replaced by a 'wand marker' to extend the marker away from the body. In this case it must be positioned to lie in the plane formed by the ASIS and PSIS

#### *Leg markers*

<b>LKNE</b>	Left knee	Placed on the lateral epicondyle of the left knee
-------------	-----------	---

Make sure the knees are as straight as possible when placing the marker (and also in the static trial). Place marker on femoral epicondyle, mid-way between anterior and posterior surfaces of the leg (excluding the patella). Palpate up from the joint line (use the popliteal crease as an indicator of the joint line, particularly in over-weight patients). Check symmetry of marker heights from the front (taking into account any leg length discrepancy). If the patient is unable to extend their knee (ie in crouch position), then place the marker slightly more posteriorly. Beware of ITB/TFL, particularly in adult males – don't allow this to influence marker positioning.

<b>LTHI</b>	Left thigh wand marker	Place the marker over the lower lateral 1/3 surface of the thigh, just below the swing of the hand, although the height is not critical.
-------------	------------------------	--

The thigh markers are used to calculate the knee flexion axis location and orientation. Place the marker over the lower lateral 1/3 surface of the thigh, just below the swing of the hand. The marker should be aligned in the plane that contains the hip and knee joint centers and the knee flexion/extension axis.

<b>LANK</b>	Left ankle	Placed on the most prominent point of the lateral malleolus along an imaginary line that passes through the transmalleolar axis. Use small base marker.
<b>LTIB</b>	Left tibial wand marker	This shank marker determines the alignment of the ankle flexion axis and is aligned by the software. Place the marker over the lower 1/3 lateral surface of shank approximately in the plane that contains the knee and ankle joint centres – i.e. lined up with the bi-malleolar axis, and not necessarily in line with LKNE.
<b>LTOE</b>	Left toe	Placed between the distal ends of the second and third metatarsals.
<b>LHEE</b>	Left heel	Placed in the middle of the calcaneus at the same height above the plantar surface of the foot as the toe marker. LTOE and LHEE used to calculate line of the foot for foot progression. <i>Note: if measuring subject in AFOs the heel marker should be in line with the heel of the foot which may be different to the heel of the AFO if there is a heel raise within the orthosis.</i>

The heel marker should be placed on the posterior aspect of the calcaneus, at the same height above the plantar surface of the foot as the toe marker. It should be aligned such that the line joining the heel and toe marker defines the progression of the foot in the transverse plane. This can be adjusted in the software if the subject is able to stand with their foot flat on the floor. If the heel is not flat on the floor this distance will need to be carefully measured.

<b>LMMA</b> (needed for static trial only)	Left medial malleolus	Placed on the most prominent point of the medial malleolus along an imaginary line that passes through the transmalleolar axis. Use small base marker.
---	-----------------------	--

Put both LANK and LMMA markers on simultaneously. This is used to define the knee/ankle axis, which is not necessarily the same as tibial torsion.

#### ***Anthropometric Measurements Required for Lower Body (to be entered on PRF form)***

<b>Height in standing (mm)</b>	In the presence of significant hip and knee contractures measure 'true' height with tape measure in supine lying).
<b>Weight (kg)</b>	taken from force plates
<b>Leg Length (mm)</b>	Distance between ASIS and the most prominent point of the medial malleolus, via the medial aspect of the knee joint. Measure in supine lying. Method takes into account any contractures and is NOT therefore the shortest distance between the ASIS and medial malleolus.
<b>Inter-ASIS distance (mm)</b>	ASIS-ASIS distance - is the distance between the left ASIS and right ASIS.
<b>Knee Width (mm)</b>	The medio-lateral width of the knee across the line of the knee axis
<b>Ankle Width (mm)</b>	The medio-lateral distance across the malleoli. Note: if measuring subject in AFOs this should be taken as the external distance including the AFOs

#### ***Optional Subject Measurements for the Lower Body***

If not entered this is calculated by the model

<b>ASIS-Trochanter Distance (mm)</b>	Vertical distance, in the sagittal plane, between the ASIS and greater trochanter when the subject is lying supine. Measure with the femur rotated such that the greater trochanter is positioned as lateral as possible. Only needed when markers cannot be placed directly on the ASIS.
--------------------------------------	---

## UPPER BODY PLUG-IN-GAIT (UL PiG)

Not routinely used in OGL clinical protocol

### Head Markers

<b>LFHD</b>	Left front head	Located approximately over the left temple
<b>RFHD</b>	Right front head	Located approximately over the right temple
<b>LBHD</b>	Left back head	Placed on the back of the head, roughly in a horizontal plane of the front head markers
<b>RBHD</b>	Right back head	Placed on the back of the head, roughly in a horizontal plane of the front head markers

The markers over the temples define the origin, and the scale of the head. The rear markers define its orientation. If they cannot be placed level with the front markers, and the head is level in the static trial, tick the "Head Level" check box under options on "Run static model" in the pipeline when processing the static trial. A headband with permanently attached markers may also be used.

### Torso Markers

<b>C7</b>	7th cervical vertebrae	Spinous process of the 7th cervical vertebrae
<b>T10</b>	10th thoracic vertebrae	Spinous Process of the 10th thoracic vertebrae
<b>CLAV</b>	Clavicle	Jugular Notch where the clavicles meet the sternum
<b>STRN</b>	Sternum	Xiphoid process of the Sternum
<b>RBAK</b>	Right back	Placed in the middle of the right scapula. This marker has no symmetrical marker on the left side. This asymmetry helps the auto-labeling routine determine right from left on the subject.

C7, T10, CLAV, STRN define a plane hence their lateral positioning is most important.

### Arm Markers

<b>LSHO</b>	Left shoulder	Placed on the Acromio-clavicular joint
<b>LUPA</b>	Left upper arm	Placed on the upper arm between the elbow and shoulder markers. Should be placed asymmetrically with RUPA
<b>LELB</b>	Left elbow	Placed on lateral epicondyle approximating elbow joint axis
<b>LFRA</b>	Left forearm	Placed on the lower arm between the wrist and elbow markers. Should be placed asymmetrically with RFRA
<b>LWRA</b>	Left wrist A	Left wrist bar thumb side, along wrist axis
<b>LWRB</b>	Left wrist B	Left wrist bar pinkie side, along wrist axis
<b>LFIN</b>	Left fingers	Actually placed on the dorsum of the hand just below the head of the second metacarpal

### Required Subject Measurements for the Upper Body

<b>Shoulder Offset (mm)</b>	vertical offset from the base of the acromion marker to shoulder joint centre
<b>Elbow Width (mm)</b>	Width of elbow along flexion axis (roughly between the distal epicondyles of the humerus)
<b>Wrist Width (mm)</b>	Anterior/ Posterior thickness of wrist at position where wrist markers are attached
<b>Hand Thickness (mm)</b>	Anterior/ Posterior thickness between the dorsum and palmar surfaces of the hand.

## OXFORD TRUNK MODEL (OTM)

The following markers are used in conjunction with the Plug-in-Gait lower body marker set.

<b>T2</b>	2 <sup>nd</sup> thoracic vertebrae	Spinous process of the 2nd thoracic vertebrae
<b>T6</b>	6th thoracic vertebrae	Spinous Process of the 6th thoracic vertebrae
<b>T12 (static only)</b>	12th thoracic vertebrae	Spinous Process of the 12th thoracic vertebrae
<b>CLAV</b>	Clavicle	Jugular Notch where the clavicles meet the sternum

## OXFORD FOOT MODEL (OFM)

The following markers are used in conjunction with the Plug-in-Gait lower body marker set.

### *Tibia markers*

<b>LHFB</b>	Left head of fibula	Most lateral aspect of the head of fibula
<b>LTUB</b>	Left tibial tuberosity	Most anterior aspect of the tibial tuberosity
<b>LSHN</b>	Left shin	Anywhere on the anterior crest of the tibia

### *Hindfoot markers*

<b>LCPG</b>	Left posterior calcaneus wand marker	Base is placed between LPCA and LHEE & its marker is in line with LPCA and LHEE, in line with coronal orientation of the hindfoot
<b>LPCA (static only)</b>	Left proximal calcaneus	Posterior aspect of calcaneus, proximal end of midline in sagittal plane
<b>LLCA</b>	Left lateral calcaneus	Lateral aspect of calcaneus at same distance from LHEE as LSTL and at same level from plantar surface of the foot as LSTL
<b>LSTL</b>	Left sustentaculum tali	Sustentaculum tali

### *Forefoot markers*

<b>LP1M</b>	Left base of 1 <sup>st</sup> metatarsal	Proximal head of 1 <sup>st</sup> metatarsal, just medial to the extensor hallucis longus tendon (palpated by asking subject to dorsiflex hallux)
<b>* LP5M</b>	Left base of 5 <sup>th</sup> metatarsal	Laterally over the base of the 5 <sup>th</sup> metatarsal
<b>* LD1M (static only)</b>	Left head of 1 <sup>st</sup> metatarsal	Medially on the head of the 1 <sup>st</sup> metatarsal
<b>* LD5M</b>	Left head of 5 <sup>th</sup> metatarsal	Laterally on the head of the 5 <sup>th</sup> metatarsal
<b>* LHLX</b>	Left base of hallux	Medial side of the proximal phalanx of the hallux, mid-way between the inferior and superior surfaces

L1DM and L5DM are positioned such that their centres lie along the line through the distal heads of the 1<sup>st</sup> and 5<sup>th</sup> metatarsal heads.


\* LP5M, LD5M, LD1M and LHLX should all be equidistant from the plantar surface of the forefoot.

Functional markers should be placed in a triangular pattern to avoid co-linearity and placed as far apart as able within the confines of available space. This would normally take the format of a distal and proximal marker on the anterior thigh and a marker at the midpoint between the knee and thigh wand. The distal thigh marker should be placed above the patella so as not to be affected by knee motion.

### **Marker cleaning**

Markers are cleaned once a week (or more regularly if required) by soaking them in hot water with washing up liquid; they are then rinsed and dried.

## Images: PiG (lower limb) with functional markers



### Functional Calibration Markers:

#### Markers to add to PiG model

R/L thigh anterior 1 = **RTHA1** & **LTHA1**

R/L thigh anterior 2 = **RTHA2** & **LTHA2**

R/L thigh lateral = **RTHL** & **LTHL**

R/L OFM fib. head = **RHFB** & **LHFB**

R/L OFM tib. tub. = **RTUB** & **LTUB**

L OFM shin = **RSHN** & **LSHN**

#### Markers to add to Oxford Foot Model

R/L thigh anterior 1 = **RTHA1** & **LTHA1**

R/L thigh anterior 2 = **RTHA2** & **LTHA2**

R/L thigh lateral = **RTHL** & **LTHL**

—————→ PiG markers

—————→ Functional markers

-----→

-----→

## Images: PiG (lower limb) with functional markers




## VST marker order


● SACR	● RTHI
● LASI	● RTHA1
● RASI	● RTHA2
● LTHI	● RTHL
● LTHA1	● RKNE
● LTHA2	● RHFB
● LTHL	● RTUB
● LKNE	● RTIB
● LHFB	● RSHN
● LTUB	
● LTIB	
● LSHN	

PiG Upper Limb (with extra UPA

and FRM) and Oxford Trunk







## Oxford Foot Model (OFM)



MARKERS IN RED ARE REMOVED AFTER STATIC TRIAL